

SCALEX[®] 800

ULTRASONIC SCALER



SCALEX[®] 800

- 25,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim lightweight handpiece.
- Convenient foot control.
- No power loss at maximum pressure.
- Unit is available in 115V 60Hz, or 230V 50Hz operation.
- Dimensions W220 x D190 x H80 mm, weight 3.2kg.


ISO 13485 **GMP**

DENTAMERICA[®]

SCALEX® 800

ULTRASONIC SCALER

✧ ULTRASONIC INSERTS



151MT(F4)

Universal Insert 25K
(External Water)

151FT(F1)

Universal Insert 25K
(Internal Water)

152MT(F5)

Chisel Type Insert 25K
(External Water)

152FT(F2)

Chisel Type Insert 25K
(Internal Water)

151CT

Endodontic Insert 25K
(Canalsonic Tip)

SPECIFICATIONS :

- * INPUT : 115V 60Hz, or 230V 50Hz
- * OUTPUT : AC 24V
- * FREQUENCY : 25,000 scaling stroke per second
- * DIMENSION : W220 × D190 × H80mm
- * WEIGHT : 3.2kg



Foot Control

DENTAMERICA®

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.asia

DISTRIBUTOR