

# SCALEX<sup>®</sup> 830

## ULTRASONIC SCALER



SCALEX<sup>®</sup> 830

- 30,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation.
- No power loss at maximum pressure.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit for 25K scaling strokes is available.
- Unit is available in 115V 60Hz or 230V 50Hz operation.

ISO 13485 **GMP**

**DENTAMERICA<sup>®</sup>**

# SCALEX® 830

## ULTRASONIC SCALER

### ※ULTRASONIC INSERTS

#### 151MT (F4)

Universal Insert 25K  
(External Water)




#### 153MT (F6)

Universal Insert 30K  
(External Water)




#### 151FT (F1)

Universal Insert 25K  
(Internal Water)



#### 153FT (F3)

Universal Insert 30K  
(Internal Water)




#### 152MT (F5)

Chisel Type Insert 25K  
(External Water)




#### 154MT (F7)

Chisel Type Insert 30K  
(External Water)




#### 152FT (F2)

Chisel Type Insert 25K  
(Internal Water)



#### 154FT (F8)

Chisel Type Insert 30K  
(Internal Water)



#### 151CT

Endodontic Insert 25K  
(Canalsonic Tip)



#### Reorder NO. A156-26

Foot control



Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

#### SPECIFICATIONS :

- INPUT : 115V 60Hz or 230V 50Hz
- OUTPUT : AC 24V 3A
- FREQUENCY : 25,000 or 30,000 Hz/sec.
- DIMENSION : W140xD140xH50mm
- WEIGHT : 1000g

### DENTAMERICA®

18688 E.San Jose Ave.  
City of industry, CA 91748 U.S.A.  
Phone:(626)912-1388  
Fax:(626)913-0510  
www.dentamerica.com  
www.dentamerica.asia

DISTRIBUTOR