

# MICROMAX® *BRITE*

Cordless Prophylaxis  
Handpiece with LED Illumination



- Built-in LED illuminates the oral cavity for easy operation.
- Cordless convenience offers complete portability.
- Easy to read LCD display.
- Heavy-duty high torque motor.
- Automatic forward/reverse functions.
- Switch from rotary to reciprocating with a press of a button.
- Nine rotation speeds available for different applications.
- High capacity lithium-ion battery.
- Four different attachments available choose from.

# MICROMAX<sup>®</sup> BRITE

Cordless Prophylaxis Handpiece with LED Illumination



## Optional Attachments



Latch Attachment



Screw Attachment



Snap-on Attachment



Straight Attachment

## SPECIFICATIONS

INPUT	AC 100V~240V, 50/60Hz
FUSE	1A 250V
BATTERY TYPE	3.6V 2040mAh
DIMENSIONS (HANDPIECE)	L201 xW26 xH31 mm
DIMENSIONS (CHARGING STATION )	L105 xW85 xH59 mm
WEIGHT (HANDPIECE)	160g
WEIGHT (CHARGING STATION )	125g

## DENTAMERICA<sup>®</sup>

18688 E. San Jose Ave.  
City of Industry, CA 91748 U.S.A.  
Phone: (626) 912-1388  
Fax: (626) 913-0510  
www.dentamerica.com  
www.dentamerica.asia

DISTRIBUTOR