



ROTEX[®] 782

COMPACT ELECTRIC HANDPIECE UNIT

ROTEX[®] 782

- Ball bearing high torque electric motor.
- 1,000 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Standard package includes Variable Speed Foot Control, Handpiece and Handpiece Rest.
- Circuit overload protection.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.
- Unit is available for 100/110V 50/60Hz or 220/240V 50Hz operation.

DENTAMERICA[®]

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.asia